

Mouse Embryonic Stem Cell Differentiation to Hematopoietic Precursors

Hogune Im*

Molecular and Cellular Pharmacology Program, Department of Pharmacology, University of Wisconsin Medical School, Madison, USA

*For correspondence: hoguneim@stanford.edu

[Abstract] Embryonic stem cells are derived from inner cell mass of an embryo that can differentiate into every cell type in the body. Clinically, cultured red blood cell supply is of great interest. However, some of the hurdles need to be overcome. This protocol describes a two step protocol to form embryoid body and differentiate them into hematopoietic lineage cells, especially erythroid cells.

Materials and Reagents

1. Mouse embryonic stem cell (ES cells)
2. 200 mM L-Glutamine (100 ml) (Life Technologies, Gibco[®], catalog number: 320-5030AG)
3. Methylcellulose (STEMCELL Technologies, catalog number: 03134)
4. Gelatin (STEMCELL Technologies, catalog number: 07903)
5. Transferrin (Human) (Boehringer-Mannheim) (Roche Diagnostics, catalog number: 652202)
6. MTG (Sigma-Aldrich, catalog number: M6145)
7. Plasma derived serum (PDS) (Animal technologies) special order from company. Typically serum contains some PDGF, which inhibit differentiation (Must call company, animal technologies, Inc)
8. Epo (Amgen) (Thermo Fisher Scientific, catalog number: 50948385)
9. Protein Free Hybridoma Media II (PFHM-II) (Life Technologies, Gibco[®], catalog number: 12040-077)
10. Ascorbic acid
11. Penicillin/Streptomycin (Pen/Strep) (Thermo Fisher Scientific, catalog number: SV30010)
12. Fetal bovine serum (FBS)
13. Fetal calf serum (FCS)
14. Phosphate buffered saline (PBS)
15. Trypsin
16. EDTA
17. NaHCO₃

18. Dulbecco's modified eagle medium (DMEM)
19. ES-IMDM (see Recipes)
20. Differentiation medium (see Recipes)
21. Secondary differentiation media (see Recipes)
22. Iscove's modified dulbecco medium (IMDM) (Life Technologies, Gibco[®], catalog number: 12200-036, 1L/pack) (see Recipes)
23. Cellulase (Sigma-Aldrich, catalog number: C-1794) (see Recipes)

Equipment

1. Standard table-top centrifuges
2. T25 flask
3. Snap cap tube
4. Incubator
5. Water bath
6. Gauge needle
7. 50 ml Falcon tube
8. 0.22 μ m filter

Procedure

- A. Primary differentiation of ES cells (making Embryoid Bodies)
 1. Two days prior to setting up differentiation, split cells into ES-IMDM medium without feeder cells in gelatinized T25 flask (this dilutes out most of the STO cells). Add 4×10^5 ES cells per T25 flask.
 2. Change medium the next day (ES cells grow faster therefore use up nutrients faster; doubles every 8 h).
 3. Set up differentiation as follows.
 - a. Aspirate medium from the flask.
 - b. Add 1 ml of trypsin/EDTA, swirl around, and quickly remove.
 - c. Add 1 ml of trypsin/EDTA and wait until cells start to detach. It usually takes about 1-2 min. Do not over-trypsinize cells.
 - d. Deactivate the trypsin by adding 1 ml FCS (differentiation serum) and 4 ml IMDM and pipette up and down to make single cell suspensions. It is important not to have cell clumps. Transfer to 15 ml snap cap tube.
 - e. Centrifuge for 5-10 min at 1,000 rpm.

- f. Resuspend the pellet in 10 ml IMDM (w/o FCS; this removes LIF and any other factors that might interfere with differentiation). Spin down cells for 5-10 min at 1,000 rpm.
 - g. Resuspend the pellet in 5 ml IMDM (w/ differentiation serum) and count the VIABLE ES cell number using eosin (0.2% eosin in 1x PBS). Make sure to count ES cells only.
 - h. Add 5,000-6,000 ES cells/ml of differentiation media for day 2.75-3 EBs. For ChIP purpose, I prepared 8-10 dishes and prepared enough media for one extra dish. Add 4,000-5,000 cells/ml for day 4-5 EBs. Add 500-2,000 cells/ml for day 6-10 EBs. Add more cells (7-8,000 cells/ml for day 2.75-3 EBs) if they differentiate poorly. Too many cells could result in poor differentiation of EBs. You don't want the cells to use up the nutrients before they reach the certain days.
 - i. Primary differentiation to derive EBs is typically done in 10 ml/10 cm dish. Use bacterial petri dish (Valmark (good), Fisher (poor) VWR (not recommended)). Do not use tissue culture dishes because all the cells will adhere and will not form EBs. If you have a way to constantly shake your dish in the incubator this improves efficiency, otherwise shake your dish 2-3 times a day until the EBs are ready.
- B. Hematopoietic precursor analysis (secondary differentiation)
1. Embryoid body (EB) harvest.
 - a. EBs in liquid: Transfer media containing EBs into 50 ml tubes. Wash the dish with IMDM and pool. Let it sit at room temperature (RT) for about 10-20 min. EBs will settle down to the bottom of the tube.
 - b. EBs in methylcellulose: Add equal volume of cellulase (2 units/ml, final 1 unit/ml) and incubate 20 min at 37 °C. Collect EBs in 50 ml tubes. Wash the dish with IMDM. Let it sit at RT for about 5-10 min. EBs will settle down to the bottom of the tube.
 2. While the EBs are settling down prepare replating media enough for one extra plate (if you need 10 plates prepare media for 11 plates) look at the recipe for details.
 3. Aspirate off as much media as possible minimizing the EB loss, add 3 ml trypsin and incubate for 3 min at 37 °C (use water bath). Very important to keep the time. It is directly related with survival and differentiation potential of ES cells. If you are harvesting for FACS analysis, 7.5 mM EDTA should be used to dissociate the EBs. This is more harsh treatment, however, receptors are kept intact.
 4. Quickly vortex briefly and add 1 ml FCS (differentiation serum).
 5. Dissociate by passing 8-10 times through a 20 gauge needle. Transfer to 15 ml tube and spin for 5-10 min at 1,000 rpm. Use collagenase for older EBs (*i.e.* >day 9). When collagenase is used, incubate cells for 1 h at 37 °C.

6. Resuspend the pellet in 0.3-1 ml of IMDM (w/ 10% FCS). Smaller the volume the better. Count the viable cells. At this point, there should be no cell clumps.
7. Add EBs to methylcellulose replating media. Shake well and let it sit for few minutes for bubbles to rise up. Using the 16 gauge blunt-end needle dispense the media+EBs on to appropriate dishes (you can use any bacterial dish for this). Always prepare bit more than what you need. Methylcellulose is sticky.
 - a. Small scale differentiation: Add $3-6 \times 10^4$ EBs/ml of methylcellulose replating media. Add 1.25 ml methylcellulose mix to each 35 mm bacterial dishes for blast colony assay. Otherwise, add 1 ml each. Prepare three replica dishes for each sample, *i.e.* make 4.5 ml methylcellulose replating media for blast colony assay and make 4 ml methylcellulose replating media for erythroid and other myeloid colony assay.
 - b. For ChIP: Add $1-1.5 \times 10^5$ EBs/ml of methylcellulose replating media. Add 10 ml methylcellulose mix to each 10 cm dishes. 10 x 10 cm plates will yield enough cells for 4 IPs ($8-10 \times 10^7$).
8. Place the plates in 245 mm square dish and a small 35 mm dish with water to maintain the humidity. This prevents the gel from drying up and cracking.

C. Harvest

1. Prepare cellulase + PBS solution (2 µg/ml)
2. Add same volume of cellulase solution to each plate so the final concentration of cellulase is 1 µg/ml.
3. Incubate at 37 °C for 30-40 min until methocell is watery.
4. Scrape plates and collect cells in 50 ml Falcon tube. Vortex and spin at 1,000 rpm, 5 min.
5. Wash with fresh media or PBS and spin at 1,000 rpm, 5 min.

Recipes

1. IMDM
 - For 1 L
 - 1 pack IMDM powder
 - 10 ml/L Pen/Strep
 - 3.024 g/L NaHCO₃
 - Filter thru 0.22 µm filter
 - If using premade liquid media -> add 5.05 ml of P/S
2. ES-IMDM

This media is essentially the same as ES-DMEM but use IMDM (richer media) instead of DMEM. Used for gelatin+ES cell culture. This is used for gelatin+ES cells and for primary differentiation to generate EBs.

3. Differentiation medium

15% pre-selected FCS in IMDM with the followings

1% ascorbic acid (5 mg/ml): Make ascorbic acid fresh each time you set up a differentiation. Dissolve ascorbic acid (5 mg/ml) in H₂O and filter sterilize (0.22 μm).

1% L-glutamine: Once thawed use it for 1 wk.

3 μl/ml MTG (MTG is diluted by adding 26 μl of MTG to 2 ml IMDM).

4. Cellulase

2 U/ml in 1× PBS filter sterilize with 0.45 μm filter.

5. Secondary differentiation media (methylcellulose mix)^a (all % is v/v)

	Erythroid	115 ml for 10 plates
Methylcellulose ^b	44%	46 ml
PDS	10%	11.5 ml
Ascorbic acid (5 mg/ml)	0.25%	287 μl
L-glutamine	1%	1.15 ml
Transferrin	1%	1.15 ml
MTG (26 μl per 2 ml) ^c	3 μl/ml	345 μl
Epo	2 unit/ml	5.75 μl
PFHM-II ^d	5%	5.75 ml
IMDM	To 100%	48,8 ml

- If efficiency of EB plating decreases, replace one or more of the following with fresh aliquots or lots: PFHM II, IMDM, MTG, FBS.
- Prewarm methylcellulose to 37 °C to decrease viscosity. Discard after it has been thawed and stored at 4 °C for >1 month. If desired, methylcellulose can be spun down at 6,000 rpm for 10 min to remove debris (transfer to new conical when spin is complete). Methocell you get from stem cell tech. is 2.5%. Calculations are made so that final methocell is 1.1%.
- Close MTG and PFHM II caps immediately after use to prevent oxidation. When opening a new 50 ml aliquots of PFHM II, sub-aliquot into five 10 ml aliquots to prevent oxidation (cover conical tubes with foil to protect from light).
- Various Factors added include bFGF (10 ng/ml), activin A (2-10 ng/ml), BMP-2 (2-10 ng/ml), and BMP-4 (2-10 ng/ml). Add KL (final 1%) and IL-3 (final 1%) for >day 6

EBs. Feed for > day 9 EBs on day 6 (4-5 ml per 10 cm petridish). Make the same methyl cellulose cocktail with 20% methylcellulose instead of 44%.

Acknowledgments

This protocol was adapted from Im *et al.* (2003). Funding for this work was obtained from the NIH.

References

1. Im, H., Park, C., Feng, Q., Johnson, K. D., Kiekhäfer, C. M., Choi, K., Zhang, Y. and Bresnick, E. H. (2003). [Dynamic regulation of histone H3 methylated at lysine 79 within a tissue-specific chromatin domain.](#) *J Biol Chem* 278(20): 18346-18352.
2. Zhang, W. J., Chung, Y. S., Eades, B. and Choi, K. (2003). [Gene targeting strategies for the isolation of hematopoietic and endothelial precursors from differentiated ES cells.](#) *Methods Enzymol* 365: 186-202.